

DLA for amplifying *Mu* flanking sequence

Digestion-ligation-amplification (DLA), a novel adaptor-mediated PCR-based method that uses a single-stranded oligo as the adaptor, was developed to overcome difficulties of amplifying unknown sequences flanking known DNA sequences in large genomes (Liu *et al.* 2009). If you have questions or comments about this protocol please contact Dr. Schnable at schnable@iastate.edu. The citation for this protocol is Liu et al., 2009 (see below).

1. Digest genomic DNA from maize plants that carry a *Mu*-induced mutant with *Nsp I* (R*CATG^Y)

1,000 ng gDNA
4 μ l 10x NEB#2 Buffer
0.4 μ l BSA (100mg/ml)
1 μ l *Nsp I* (NEB, Cat#:R0602L, 10,000U/ml) = 10U
ddH₂O to 40 μ l

Incubate at 37°C for 1.5hr

2. Ligation using adaptor

Prepare 50uM *NspI*-5B Adaptor Mix (100ul):

| | |
|-------------------------------|------|
| 100uM <i>NspI</i> -5B*(oligo) | 50ul |
| 10X NEB#2 buffer | 5ul |
| ddH ₂ O | 45ul |

* see below for primer sequences

After 1.5hrs incubation, add to genomic DNA digestion reaction:

5.0 μ l *NspI*-5B Adaptor Mix (50uM)
3.0 μ l ddH₂O
1.0 μ l 10X Ligase Buffer
1.0 μ l T4 ligase (NEB, Cat# M0202L)

Total volume per reaction is now 50 μ l. Incubate for 3 hrs at 16°C in the thermocycler (Program "Dig-Lig", see below).

3. Purification – Follow the protocol provided for PCR purification with the Qiaquick PCR purification kit (QIAGEN, Cat# 28106 or Cat# 28104).

4. (Optional) Blocking using ddNTP

| | |
|--------------------------|------------|
| Purified Dig-Lig product | 50-200ng |
| Klenow 10X buffer | 5 |
| 2mM ddNTP | 2 |
| Klenow 5U/ul | 0.8 |
| ddH ₂ O | Up to 50ul |

(See below for the program: Klenow30)

Note: This is an optional step. You can try to use the protocol without this step first. If too much background amplification occurs, repeat with the addition of this step.

5. Purification – Follow the protocol provided for PCR purification with the Qiaquick PCR purification kit (QIAGEN, Cat# 28106 or Cat# 28104).

CAUTION: It's important to fully remove ddNTP. It is therefore suggested to conduct the PE wash step twice.

6. 1st PCR (PCR program: HL60-30, see below.)

| | |
|----------------------------|----------|
| GeneAmp 10X PCR Buffer II | 2ul |
| 2mM dNTP | 2ul |
| MgCl ₂ (25mM) | 1.2ul |
| MuTIR primer (5uM)* | 1.2ul |
| NspI-5 primer (5uM)* | 1.2ul |
| Purified DNA (from step 3) | 50~100ng |
| AmpliTaq Gold (5U/ul)** | 0.2ul |
| Water | To 20ul |

* see below for primer sequences

** AmpliTaq Gold® DNA Polymerase (Applied Biosystems, Cat#: N808-0241)

7. Nested PCR (PCR program: HL60-35, see below)

| | |
|---|-------|
| GeneAmp 10X PCR Buffer II | 3ul |
| 2mM dNTP | 4ul |
| MgCl ₂ (25mM) | 1.8ul |
| Mu53s* (5uM) | 1.5ul |
| NspI-P* (5uM) | 1.5ul |
| 10X diluted 1st PCR product (from step 4) | 2ul |

| | |
|-----------------------|---------|
| AmpliTaq Gold (5U/ul) | 0.2ul |
| Water | To 30ul |

* see below for primer sequences

8. Oligo sequences:

- 1) MuTIR 32 mer
5' AGAGAAGCCAACGCCAWCGCCTCYATTTTCGTC
- 2) Mu53s 19 mer
5' GCCTCYATTTTCGTCGAATC
- 3) NspI-5B 22 mer
5' CAGAACGTCACAGCATGTCATG
- 4) NspI-5 20 mer
5' GAACGTCACAGCATGTCATG
- 5) NspI-P 19 mer
5' AACGTCACAGCATGTCATG

9. PCR programs

- 1) Dig-Lig
1.5h @ 37°C
4°C
3h @ 16°C
4°C forever
- 2) Klenow30
0.5h @ 30°C
75°C 10min
4°C forever
- 3) HL60-30
94°C 10min
30 cycles { 94°C 30s
60°C 45s
72°C 2.5min
72°C 10min
4°C forever
- 4) HL60-35
94°C 10min
35 cycles { 94°C 30s
60°C 45s

72°C 2.5min
72°C 10min
4°C forever

References

LIU, S., C. R. DIETRICH and P. S. SCHNABLE, 2009 DLA-Based Strategies for Cloning Insertion Mutants: Cloning the *gl4* Locus of Maize using *Mu* Transposon Tagged Alleles. Genetics **in press**.